
Aesculap Power Systems


-
- (GB)** **Instructions for use/Technical description**
macro-Line transversal saw GB124R
 - (D)** **Gebrauchsanweisung/Technische Beschreibung**
macro-Line Transversalsäge GB124R
 - (F)** **Mode d'emploi/Description technique**
Scie transversale macro-Line GB124R
 - (E)** **Instrucciones de manejo/Descripción técnica**
Sierra transversal macro-Line GB124R
 - (I)** **Istruzioni per l'uso/Descrizione tecnica**
Sega trasversale macro-Line GB124R



GC523


GC500

GC502

GC510

GA550—GC554

macro-Line transversal saw GB124R

Legend

- 1 Double wrench
- 2 Nut
- 3 Lock washer
- 4 Disk (with two slits)
- 5 Saw blade slot
- 6 Saw blade
- 7 Pin
- 8 Triangular shaft
- 9 Denticulation
- 10 Lubrication point (for clamping lever)
- 11 Fixing screw
- 12 Clamping lever
- 13 Chuck
- 14 Tool coupling
- 15 Fixing screw
- 16 Hexagon wrench
- 17 Face plate (motor unit)
- 18 Nose
- 19 Centering lug
- 20 Triangular shaft

Symbols on product and packages

Symbol	Explanation
	Caution: See documentation supplied with the product

Contents

1. Safe handling	3
2. Product description	3
2.1 System components	3
2.2 Components necessary for use	3
2.3 Intended use	3
2.4 Operating principle	3
3. Preparation and setup	3
4. Working with the macro-Line transversal saw GB124R	3
4.1 System set-up	3
Connecting the accessories.....	3
4.2 Function checks	3
Test run	3
4.3 Safe operation	4
Coupling a motor unit	4
Uncoupling the motor unit	4
Mounting a saw blade on the shaft	4
Clamping a saw blade.....	4
Sawing	4
5. Validated processing procedure	5
5.1 General notes	5
5.2 Preparations at the place of use	6
5.3 Preparation prior to cleaning	6
5.4 Cleaning/Disinfecting	6
5.5 Manual cleaning/disinfecting	7
Manual cleaning and wipe disinfection	7
5.6 Inspection, maintenance and checks	7
5.7 Packaging	7
5.8 Sterilization method and parameters	8
5.9 Storage	8
6. Maintenance	8
7. Troubleshooting list	8
8. Technical Service	9
9. Accessories/Spare parts	9
10. Technical specifications	9

1. Safe handling

- Thoroughly clean the new product after removing its transport packaging and prior to its initial sterilization.
- Prior to use, check for proper condition and functioning of the product.
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Observe all safety information and maintenance advisories.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Keep the instructions for use accessible for the user.

2. Product description

2.1 System components

Designation	Art. no.
Handpiece	GB124R
Hexagon wrench	GB020R
Double wrench (SW 14)	GC523

2.2 Components necessary for use

- macro-Line transversal saw GB124R
- Flexible shaft GA156 (1.6 m) with Elan® EC and foot switch
 - or –
- Flexible shaft GA172 (2.1 m) with Elan® EC and foot switch
 - or –
- AIRLAN® pneumatic motor GA200

2.3 Intended use

The macro-Line transversal saw GB124R is used for the dissection of bones in osteosynthesis surgery.

2.4 Operating principle

The saw handpiece provides 1:1 conversion of the rotational movement of the motor unit to oscillating movement.

The bones are dissected by the oscillating movement of the saw blade. Any soft tissue parts touched by the saw blade remain unharmed.

3. Preparation and setup

Non-compliance with these rules will result in complete exclusion of any responsibility and liability on the part of Aesculap.

4. Working with the macro-Line transversal saw GB124R

4.1 System set-up

Connecting the accessories

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

4.2 Function checks

- Check all cutting edges for any mechanical damage.

Test run

- Briefly run the transversal saw at maximum speed. Check for proper and reliable functioning.

macro-Line transversal saw GB124R

4.3 Safe operation

WARNING

Risk of burns to patient and user, due to hot motor/hot tool!

- Put down the motor/hot tool beyond reach of the patient.
- Use a cloth to protect against burns when changing the motor or tool.

WARNING

Risk of injury and material damage due to inappropriate use of handpieces or motors!

- Always follow the safety advice and information given in the instructions for use of the handpieces and motors.
- Do not exceed the specified maximum motor speed.

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check before using the product.

WARNING

Risk of injury when applying the product outside the field of view!

- Apply the product only under visual control.

Coupling a motor unit

The saw can be powered via a flexible shaft GA156/GA172 or by the AIRLAN® pneumatic motor GA200.

- Turn clamping lever **12** towards the **AUF** (OPEN) position and release the clamping of the transversal saw.
- Install the handpiece on centering lug **19** in such a way that it touches the face plate **17** of the motor unit and nose **18** engages in the recess at the handpiece.
- Turn clamping lever **12** towards the **ZU** (CLOSED) position.

Moving the clamping lever to be parallel with the motor unit

- Tighten clamping lever **12** hand-tight.
- Undo fixing screw **11** with hexagon wrench **16** until the denticulation **9** disengages.
- Move clamping lever **12** to the intended position and engage denticulation **9**.
- Tighten fixing screw **11** with hexagon wrench **16**.

Uncoupling the motor unit

- Turn clamping lever **12** towards the **AUF** (OPEN) position and release the clamping of the transversal saw.
- Remove the motor unit.

Mounting a saw blade on the shaft

- Install saw blade **6** in the intended position on shaft **8**. Be careful to position pin **7** in saw blade slot **5**.
- Install disk **4** in such a way that a slit engages over pin **7**.
- Position lock washer **3** and nut **2** and tighten the nut with double wrench **1**.

Clamping a saw blade

Note

Any saw blade or trepan with triangular shaft **20** can be clamped in.

- Insert triangular shaft **8** in tool coupling **14** and secure it with fixing screw **15**.

Sawing

The optimum cutting performance is achieved by sawing with high motor speed and moving the handpiece back and forth.

Note

To avoid resonant oscillations of the saw blade, do not run the handpiece at maximum speed without load.

- Start the sawing procedure at moderate motor speed, which is then slowly increased.

5. Validated processing procedure

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, the manufacturer's specifications can only serve as an approximate guide for assessing the processing procedures applied by the individual operator/processors.

5.1 General notes

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. To avoid this, the time interval between application and processing should not exceed 6 h, and neither fixating pre-cleaning temperatures >45 °C nor any fixating disinfecting agents (active ingredients: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only process chemicals that have been tested and approved (e.g. VAH/ DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Optical changes to the material, e.g. fading or discoloration of titanium or aluminum. For aluminum, pH >8 in the application/process solution is sufficient to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use process chemicals that cause stress cracking or brittleness of plastic materials.
- For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see www.a-k-i.org

macro-Line transversal saw GB124R

5.2 Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.

5.3 Preparation prior to cleaning

- Have the product dry in a disposal container and ready for immediate cleaning and disinfecting within 30 min after use.

5.4 Cleaning/Disinfecting

CAUTION

Damage to, or destruction of the product caused by mechanical cleaning/disinfecting!

- Only clean and disinfect the product manually.

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents!

- Only use cleaning/disinfecting agents approved for surface cleaning. Follow the manufacturer's instructions for the respective cleaning/disinfecting agent.

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to manufacturers' instructions. The cleaning/disinfecting agents must
 - be approved for plastic materials and high-grade steel,
 - not attack softeners (e.g. silicone).
 - Observe specifications regarding concentration, temperature and exposure time.
- Do not clean motors/handpieces by ultrasonic treatment and do not immerse them in any fluid.
Let any fluid that penetrates into the product drain out to avoid the risk of corrosion/malfunctioning of the product.

5.5 Manual cleaning/disinfecting

Manual cleaning and wipe disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	-	-	D-W	-
II	Drying	RT	-	-	-	-
III	Wipe disinfection	-	≥1	-	-	Meliseptol HBV wipes 50 % propan-1-ol
IV	Final rinse	RT (cold)	0.5	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Clean the product under running tap water, using a suitable cleaning brush if necessary, until all visible residues have been removed from the surfaces.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Make certain the product is positioned in such a way that water will not enter the product e.g. through coupling interfaces. (Immediately remove any fluid that entered the product inadvertently.)
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

Stage III

- Wipe all surfaces of the product with a single-use disinfecting wipe.

Stage IV

- After the specified exposure time (of at least 1 minute), rinse the disinfected surfaces under running FD water.
- Hold the product in such a way that water cannot enter it e.g. through couplings or adapters.
- Allow water to drip off for a sufficient length of time.

Stage V

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

5.6 Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functional and undamaged.
- After every cleaning and disinfecting cycle, lubricate the product at lubrications point **10** and centering lug **19** with one drop from STERILIT® Power Systems drip lubricator GA059.
- Check the product for any damage, atypical running noise, overheating or excessive vibration.
- Inspect tools for broken, damaged or blunt edges.
- Set aside the product if it is damaged.

5.7 Packaging

- Always observe the instructions for use of packaging and storage devices.
- Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product (DIN EN ISO 11607).

macro-Line transversal saw GB124R

5.8 Sterilization method and parameters

- Prior to sterilization, remove any tools and adapters and uncouple the handpiece from the motor unit.
- Check to make certain that the sterilizing agent will be in contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/2 bar; holding time 5 min.
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

5.9 Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

6. Maintenance

To ensure reliable operation, Aesculap recommends having the product serviced at least once a year.

For services to that end, please contact your national B. Braun/Aesculap agency, see Technical Service.

7. Troubleshooting list

Malfunction	Cause	Finding	Remedy
Excessive noise	Motor unit not coupled correctly	Motor unit can be pulled out from handpiece	Couple the motor unit, see Coupling a motor unit
	Worn gearing	Handpiece hot	Have product repaired by the manufacturer
Lack of power at handpiece	Worn gearing	Excessive noise	Have product repaired by the manufacturer
Insufficient cutting power of saw blades	Saw blade blunt	Saw blade worn	Replace saw blade
	–	Motor unit not running correctly	Observe the instructions for use of the motor unit
Saw blade not moving	Coupling defective	Coupling parts broken	Have product repaired by the manufacturer
	Defective gearing	Motor unit running, but not powering the handpiece	Have product repaired by the manufacturer

8. Technical Service

Risk of injury and/or malfunction!

➤ **Do not modify the product.**

WARNING

➤ For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1601
Fax: +49 7461 14-939
E-mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

9. Accessories/Spare parts

Art. no.	Designation
GA059	STERILIT® Power Systems drip lubricator
GA156	Motor unit: Flexible shaft
GA172	
GA200	AIRLAN® pneumatic motor
GB020R	Hexagon wrench
GC500R	Saw blade only GC501 and shaft GC520
GC502R	Saw blade only GC503 and shaft GC520
GC504R	Saw blade only GC505 and shaft GC520
GC510R	Saw blade only GC511, GC512 and shaft GC521
GC513R	Saw blade only GC514, limiting disk GC515 and shaft GC521
GC520/204	Nut
GC523R	Double wrench (SW 14)
GC550R	Hollow reamer, Ø: 8 mm
GC554R	Hollow reamer, Ø: 16 mm
GC562R	Trepan, Ø: 8 mm
GP024/210	Disk (with two slits)
GP024/211	Lock washer
TA004103	Fixing screw

10. Technical specifications

Classification acc. to Directive 93/42/EEC

Weight	300 g
Maximum motor speed	20 000 1/min
Transmission Rotation to oscillation	1:1
Saw blade coupling	Triangular shaft
suitable for sterilization	up to 143 °C (3 bar)

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Aesculap AG