**DE**

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## FOKUSEINSTELLUNG

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## FOCUSING SYSTEM

## REPLACING THE BATTERY

The Diascleral Transilluminator is supplied with 2 AAA batteries (1.5 V). To replace the battery, unscrew the cap at the end of the light. Since the sockets and bulbs vary in measurements due to temperatures, a grub screw is placed in the socket, fixing the bulb.

Before removing the bulb, please make sure, the grub screw is loosened by about 1 mm. The appropriate allen key is included with every Diascleral Transilluminator.

Please remove the discharged battery and dispose of it properly. Insert the new battery **with the positive pole (+) first**. Close the Transilluminator by screwing the end cap back in place. The end cap must be fully secured to ensure complete electrical connection is established.

If you do not intend to use the Transilluminator for a long period of time, remove battery to prevent it leaking and damaging the lamp.

## CLEANING

The cone is made of acrylic glass. Alcoholic detergent and disinfectants can crack and/or blur the glass. Please clean the cone with a mild detergent only.

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#### UNITED KINGDOM

EICKEMEYER Ltd.  
3 Windmill Business Village  
Brooklands Close  
Sunbury-on-Thames  
Surrey, TW16 7DY  
T +44 20 8891 2007  
F +44 20 8891 2686  
E info@eickemeyer.co.uk  
www.eickemeyer.co.uk

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