



LifeVet CP Capnograph
and Pulse Oximeter

veterinary technology for life
Eickemeyer®

Smart Technology made easy



LifeVet CP Capnograph and Pulse Oximeter

The LifeVet CP handheld Capnograph and Pulse Oximeter measures and displays the expiratory CO₂ concentration as well as the respiratory rate, pulse rate and oxygen saturation in dogs, cats and horses. This compact and portable device is very easy to use with its touchscreen display and intuitive interface.

Features

- ▶ High resolution 4.3" Touch Screen Display
- ▶ Side-stream Capnography
- ▶ EtCO₂ concentration: 0 – 150 mmHg
- ▶ Respiratory Rate displayed in rpm
- ▶ SpO₂ range: 1 – 100 %
- ▶ Pulse Rate: 25 – 250 bpm
- ▶ Audible and visual alarms with adjustable volume
- ▶ Long-life Lithium battery: Up to 18 h operating time
- ▶ Dimensions (in cm): W 7,7 x H 14,1 x D 3,6
- ▶ Weight: 420 g
- ▶ Includes: CO₂ tube, CO₂ filter, net adapter, USB cable, clip sensor
- ▶ Please note: Please order the exhaust air connection set (321881)

321880



EICKEMEYER®

Via G. Verdi 8
65015 Montesilvano (PE)
Italia

T 085 935 4078
F 085 935 9471
info@eickemeyer.it

TELEFONO 085 935 4078

www.eickemeyer.it