

DENTO-DORM® Tracheal Sponge Set

Item no. 178260

DENTO-DORM® Tracheal Sponge Made from cellular rubber, this is a sealing sponge for the trachea, especially during calculus removal and dental treatments. The DENTO-DORM® Tracheal Sponge has a tightly sealing outlet opening for the tracheal tube. Available as a basic case kit with various tracheal sponge sizes for small cats to large dogs. The case comes with 48 tracheal sponges.

Sponge seals with central, densely sealing perforations for the tracheal tube. Employed during calculus removal and dental treatments far in the pharynx and prevents calculus fragments and water from reaching the pharynx. At the same time ensures the tube from slipping.

Operating at a higher level


The idea was born during the daily routine in surgery. On the one hand, it was meant to offer the animal optimal heat-insulating lying and sleep comfort. On the other, the treating veterinarian's work was supposed to be made easier and more ergonomic, improving the results at the same time. This resulted in the completely newly developed, high quality and long-lasting DENTO-DORM® products.

Secured Quality

DENTO-DORM® is developed, produced and clinically tested in Germany according to strict standards. Thanks to physically interlaced and stable foam material with a sealed surface, we can offer our customers the highest quality.

Sophisticated Design

Our heat-insulating positioning mattresses offer highest flexibility for dental treatments, X-rays, sonographies, blood collection and surgeries on dogs, cats and small pets. They can be combined in all variations.



application