




Anaesthetic Machine
LAPD-1000

veterinary technology for life
Eickemeyer®

Sweet dreams



Anaesthetic Machine LAPD-1000

The LAPD-1000 is a passive anaesthesia machine that is the ultimate blend of simplicity, ergonomics, and efficiency.

The machine was designed to be user-friendly and to allow the operator to expand the capabilities of their practice and use it together with a ventilator. Suitable for large animals.

Includes

- ▶ Flowmeter (0 to 10 litres/min)
- ▶ TEC 3 Isoflourane Vaporiser
- ▶ 30 litre Breathing bag
- ▶ Y-Piece Connector and 4 adaptors (18/20/25/30 mm)

213101

Accessories (not shown):

Ventilator Machine LAVC-3000

The LAV-3000 is designed for facilities that currently use a passive anaesthesia machine and wish to expand their capabilities to include full ventilation of their patients. Suitable for large animal anaesthesia machines.

213120